



GUIDELINE for the forensic examination of paint by Raman spectroscopy

DOCUMENT TYPE:	REF. CODE:	ISSUE NO:	ISSUE DATE:
GUIDELINE	EPG-GUIDELINE-003	001	01.05.2022

1

2 GENERAL REMARK

3 This guideline assumes prior knowledge in the forensic discipline. It is based on consensus
4 among the relevant forensic experts and reflects the accepted practices at the time of writing.
5 The requirements of the judicial systems are addressed in general terms only.

6

7

8 CONTENTS

9

	Page
1. AIMS	1
2. SCOPE	1
3. DEFINITIONS AND TERMS	2
4. TECHNICAL BACKGROUND	3
5. SAMPLE HANDLING AND PREPARATION	6
6. MEASUREMENT CONDITIONS AND PARAMETERS	6
7. INTERPRETATION OF RESULTS	9
8. REFERENCES	10
9. AMENDMENTS AGAINST PREVIOUS VERSIONS	11

10

11

12 1. AIMS

13 This guideline provides information and recommendations on the analysis of coating materials
14 by Raman spectroscopy. It is designed to be used in conjunction to the Best Practice Manual
15 for the forensic examination of paint (EPG-BPM-001).

16

17

18 2. SCOPE

19 This guideline is intended to provide understanding of the possibilities, advantages, limitations
20 and the proper use of Raman spectroscopy.

21

22

23

24 **3. DEFINITIONS AND TERMS**

25 For the purposes of this guideline, the relevant terms and definitions are given in ENFSI
 26 documents, in ILAC G19 [1], in ISO/IEC 9000 [2], ISO/IEC 17020 [3] and ISO/IEC 17025 [4]
 27 standards, and in ASTM E1610-18 [5] are applied. Specific technical terms used in this
 28 guideline include:
 29

Term	Definition
<i>Anti-stokes lines</i>	Scattered light with higher energy (lower wavelength) than the exciting laser light
<i>Confocal</i>	Optical imaging technique used to increase micrograph contrast by using a pinhole to eliminate out-of-focus light
<i>Dispersive instrument</i>	Separating different wavelengths by using a grating-based dispersive unit (monochromator)
<i>Edge filter</i>	Band-pass filter that passes frequencies within a certain range and rejects frequencies outside this range. The edge filter rejects both the excitation laser line and the anti-Stokes lines.
<i>Fluorescence</i>	Radiation emitted from the excited sample that overlays the weak Raman signals
<i>Laser</i>	Device for emitting monochromatic light through a stimulated emission process
<i>NIR</i>	Near infrared radiation (wavelength interval of 700 to 1400 nm)
<i>FT-Raman</i>	Separating different wavelengths by measuring a time dependent interferogram and transforming it by Fourier Transformation into a frequency function (spectrum)
<i>Notch filter</i>	Band-stop filter with a narrow stop band rejecting the excitation laser line and passing both the Stokes and anti-Stokes lines
<i>Polarizability</i>	Relative tendency of a charge distribution (electron cloud) of an atom or molecule to shift under the influence of external stimulation
<i>Raman scattering</i>	Inelastic scattering of photons from a source
<i>Raman shift</i>	Difference between the excitation wavelength and the wavelengths of the scattered light
<i>Raman spectroscopy</i>	Spectroscopic method to compare and identify substances by measuring the shift of scattered light after monochromatic excitation
<i>Rayleigh scattering</i>	Elastic scattering of the excitation light by a sample
<i>SERS</i>	Surface enhanced Raman scattering is a surface sensitive technique that results in the enhancement of the Raman scattering by molecules adsorbed on a metal or metal colloid substrate, the enhancement factor can be as much as $10^{14} - 10^{15}$.
<i>SERRS</i>	Surface enhanced resonance Raman scattering is a combination of the SERS technique with a laser excitation near the absorption maximum of the analytes
<i>Stokes lines</i>	Scattered light with lower energy (longer wavelengths) than the exciting laser light

30
 31
 32

4. TECHNICAL BACKGROUND


4.1 Theoretical background

Infrared and Raman spectroscopy both measure vibrational transitions in molecules, but these methods are based on different selection rules. To be active in infrared the dipole moment of the molecule must change. For a transition to be Raman active there must be a change in polarizability of the molecule. Thus, Raman and Infrared spectroscopy are complementary.

Raman spectroscopy is based on the scattering of light by matter. When monochromatic light of energy $h\nu_0$ encounters a sample, there is a high probability that it will be elastically scattered i.e. with the same frequency (Rayleigh scattering). About 1 in 10^{10} incident photons however undergoes inelastic scattering i.e. with an energy that differs from $h\nu_0$ by an amount equal to a vibrational transition $h\nu_s$:

Stokes Raman scattering occurs when a molecule absorbs an incident photon but relaxes to a higher vibrational state instead of the ground state

Anti-stokes Raman scattering occurs when a molecule in a vibrational excited state absorbs an incident photon and relaxes to the ground state instead of the initial excited state.



Stokes and anti-Stokes peaks are symmetrically positioned about the Rayleigh peak but the former is more intense. As the effects are weak and fluorescence can overwhelm the peaks, Raman spectra are usually recorded in the Stokes region only. Monochromatic excitation is performed using laser sources.

Raman spectra consist of scattered intensity plotted against Raman shift from the incident energy, usually expressed in wavenumbers (cm^{-1}).

4.2 Instrumentation

Raman spectroscopic measurement of paint particles requires the use of a spectrometer coupled with a microscope.

70 4.2.1 Raman spectrometer

71 Two types of instruments are available:

72

73 Dispersive spectrometers separate the scattered wavelengths using a grating and
74 measure them on a thermoelectrically cooled CCD detector. Resolution of these
75 instruments depends on the grating used (number of lines per mm).

76

77 FT-Raman spectrometers use an interferometer to modulate the beam similar to FTIR
78 instruments. High sensitivity, liquid nitrogen cooled Ge detectors are used.

79

80 The choice of the laser wavelength is of prime importance for Raman spectroscopy.
81 Measuring paint samples requires that several (at least three advised) different wavelength
82 lasers should be available in order to avoid fluorescence and maximize detection capabilities
83 (especially in case of mixed samples). Laser sources in the blue-green and red regions of the
84 visible range and a NIR source are complementary and enable measuring most samples.

85

blue	457 nm	diode-pumped solid-state laser
blue	473 nm	diode-pumped solid-state laser
blue	488 nm	Ar ion (outdated)
blue	488 nm	optically-pumped diode laser
green	514 nm	Ar ion (outdated)
green	514 nm	diode-pumped solid-state laser
green	532 nm	solid-state laser
red	633 nm	He-Ne laser
red	638 nm	diode laser
NIR	785 nm	diode laser
NIR	1064 nm	Nd:YAG

86

87 Short wavelengths are usually detected with higher efficiency and less noise but are likely to
88 excite more fluorescence. FT-Raman spectrometers use a 1064 nm laser.

89

90 Lasers should be tuned on and allowed to reach stability prior to operational runs (according
91 to manufacturers' recommendations, approx. 20 min).

92

93 4.2.2 Microscope

94 Any stable reflected light microscope can be adapted to the spectrometer. A range of infinity
95 corrected objectives are used ranging from 10x to 100x. The spot size of the laser beam on
96 the sample and the depth of focus will decrease as the magnification increases (e.g. 50x
97 objective, 514 nm laser gives a spot size of 6,5 μm and depth of focus of 10 μm). They will
98 increase with the laser wavelength.

99

100 Confocal microscopes have an additional aperture (confocal hole) that prevents out-of-focus
101 rays to reach the detector. Confocal mode will enable the reduction of the depth of focus to
102 less than 2 μm (514 nm laser, 100x objective), but also reduces the Raman signal. This mode
103 can be used to minimize fluorescence or to reduce interference from surrounding materials.
104 Changing the laser wavelength is generally more efficient in minimising fluorescence.

105

106

107 4.2.3 Rejection filter
108 Rejection filters are inserted in the optical path after the sample and transmit the Raman signal
109 to the detector while rejecting the intense Rayleigh line. They are matched to one specific
110 laser wavelength and have to be changed when changing the laser excitation wavelength.

111
112 Holographic notch filters can be used but have a limited lifetime.

113
114 Oxide coated edge filters reject both the Rayleigh line and the anti-Stokes scattering. As only
115 Stokes scattering is used for measuring paint samples, edge filters are preferred.

116
117 4.3 Use in paint analysis
118

119 Many pigments and fillers used in paint formulations do not produce unambiguous information
120 in the FTIR spectra due to unfavourable selection rules and the domination by the binders
121 response. As a complementary technique to FTIR, Raman spectroscopy usually provides
122 clear pigment and filler signals as opposed to weak response for the binders. It thus provides
123 the means to compare these components and identify them using reference databases.

124
125 Different paints with the same pigment/filler composition can be indistinguishable using
126 Raman spectroscopy.

127
128 The use of a microscope coupled with a Raman spectrometer has enabled the analysis of
129 compact paint particles as well as smears. It is possible to analyse very small samples (0,01
130 mm or less) quickly, reproducibly and with a high degree of sensitivity. The diameter of the
131 spot size is usually a few microns.

132
133 The method is practically non-destructive and provides additional or at least confirmatory
134 information to FTIR analyses.

135
136 4.4 Limitations
137

138 Using a Raman microscope and the most commonly used objective 50x, the laser spot size
139 at the sample surface is 2 μm in diameter (lateral resolution) and the penetration depth is
140 approximately 2.5 μm (confocal resolution). The following limitations are given:

141

amount	:	Minimum spot size 4 μm^2
Detection limit	:	Approx. 10 pg (2 μm \varnothing x 2 μm H, density 1.5 g/cm ³)

142
143 Organic coloured pigments measured with specific laser wavelength may show resonance
144 enhancement, yielding an improved detection limit.

145
146 A drawback of the method is the possible localised thermal degradation of the paint if the laser
147 source power is not carefully managed.

148
149 A major limitation is the occurrence of fluorescence emission that can partially or totally mask
150 the detection of the Raman signal. The use of multiple laser sources helps in overcoming this
151 limitation.

152
153

154 **5. SAMPLE HANDLING AND PREPARATION**

155

156 **5.1 Sample cleaning**

157

158 Any impurity on the sample surface has to be documented, removed and, if forensically
159 significant, analysed.

160

161 Depending on the sample, impurities can be removed either by scraping with a scalpel, or
162 washing with a tissue soaked with water or a non-aggressive solvent like ethanol or methanol.
163 Other organic solvents must be avoided because of possible chemical reactions with the
164 sample or other alterations (dissolution, swelling).

165

166 **5.2 Preparation techniques**

167

168 Sample preparation techniques should be the same for all samples being compared.

169

170 Raman spectra are measured in reflectance. The samples are deposited either directly on a
171 microscope slide, on a microscope slide covered with aluminium foil, or better still, on a slide
172 of aluminium. The metal provides advantages due to its heat dissipation properties.

173

174 Smears of paint on other surfaces (e.g. metals, plastics, textiles, paint chips) can be measured
175 directly on these surfaces. In this case a comparison measurement of the pure surface is
176 necessary to eliminate signals of the latter.

177

178 In case of multilayer samples thin sections may eventually be prepared to provide selective
179 analysis per layer. The thin sections should be placed on an appropriate sample holder or
180 again on an aluminium (foiled) slide.

181

182 **5.3 Special techniques: SERS, SERRS**

183

184 SERS and SERRS enhance the Raman signal and quench fluorescence.

185

186 Several procedures exist to prepare SERRS solution, e.g. citrate-reduced silver colloid in
187 aqueous solution of poly (L-lysine).

188

189 The techniques require good contact between the reagent and the sample, which is often
190 difficult to obtain because pigments are trapped in the binder and thus are not in direct contact
191 with the reagent spotted on the paint surface.

192

193 These special techniques are not yet applicable in routine analysis of paint.

194

195

196 **6 MEASUREMENT CONDITIONS AND PARAMETERS**

197

198 **6.1 Calibration and validation**

199

200 Instrument calibration and validation records shall be maintained in order to check the
201 instrument performance and long-term stability.

202

203

204 6.1.1 Frequency calibration

205 Raman shift values are more prone to error than in FTIR systems. They have to be calibrated
206 because small changes in the true wavelength of the laser sources can have a significant
207 impact on the Raman shift accuracy.

208
209 Frequency calibration with absolute frequency standards is usually done by the manufacturer.
210 A Ne, Ar or Hg lamp is placed near the sample position and the atomic emission lines are
211 used as well as the lines from the laser itself (must correspond to zero Raman shift).

212
213 The ASTM committee collected a set of 8 Raman shift standards [6] with a standard deviation
214 of less than 1 cm^{-1} . These standards can be used independently of the laser frequency,
215 provided it is constant.

216
217 Among these standards, polystyrene is particularly suitable and is already used for FTIR
218 calibration. It shows Raman shifts between 620 and 2940 cm^{-1} that are determined using 514
219 and 1064 nm lasers.

220
221 Another Raman shift standard is silicon with a sharp single band at 520 cm^{-1} .

222
223 One of these standards should be used for accuracy check when starting a session of
224 measurements. Accuracy checks should concern laser lines and gratings that will be used for
225 measurements.

226
227 6.1.2 Performance monitoring
228 The instrument response has to be controlled by checking the magnitude of the Raman signal
229 under the same analytical conditions. This must be done for each laser source available with
230 the instrument. The intensity of the silicon band at 520 cm^{-1} thus permits the detection of
231 instrument changes (e.g. misalignment of the laser beam, output power decrease of the laser
232 source).

233
234 In addition, a secondary standard similar to the sample of interest may be chosen, e.g. a
235 known pigment from the laboratory collection or polystyrene. In this way positions and relative
236 intensities of different peaks can be controlled in the spectral range of the sample interest.

237
238 6.1.3 Validation
239 Raman shift standards can be instrumental for validation and provide measurands in the form
240 of the position and relative intensities of the peaks. The chosen peaks must cover a wide
241 range of Raman shift.

242 In addition, other standards similar to the sample of interest should be chosen, e.g. a known
243 pigment from the laboratory collection. Positions and relative intensities of different peaks must
244 be controlled.

245
246 6.2 Analysis conditions

247
248 Analytical conditions have to be optimized according to instrument and sample type.
249 Therefore only generic advice can be formulated.

250
251 6.2.1 Wavelength range
252 A range of $2000 - 200\text{ cm}^{-1}$ will be sufficient for measuring most pigments. Certain blue and
253 green pigments also display spectral characteristics up to 4000 cm^{-1} .

254

255 6.2.2 Microscope objective
256 The choice of the objective will influence the size and depth of the area measured. Different
257 objectives should be tested on the samples. Generally, 50x or 100x objectives will provide the
258 best results.

259
260 6.2.3 Laser power
261 The laser power should not saturate the detector nor burn the sample. It is recommended to
262 start with a low power and increase it depending on the result obtained. A high laser power
263 may modify the molecular structure and thus the Raman spectra of heat sensitive samples
264 (i.e. graphitization bands may appear).

265
266 6.2.4 Analysis time, number of accumulations, number of measurements
267 These parameters should be chosen so that an acceptable S/N ratio is obtained. As a start,
268 the time of analysis can be short (i.e. 10 s), no accumulation, 1 measurement. If the S/N
269 ratio is good, there is no need to adjust the parameters. Normally the same parameters are
270 used for both the control and the questioned samples. Sometimes the intensity of the
271 background fluorescence may require the use of different parameters for the control and the
272 questioned samples, or even within the same sample.
273 Acquiring at least two accumulations may be useful to detect or automatically remove
274 artefacts (if available in the software), such as 'cosmic rays' or 'spikes' in the spectrum.

275
276 6.2.5 Grating
277 Some instruments allow the grating being changed (from 300 to 1800 lines per mm). The
278 choice of the grating depends on the laser wavelength. 300 lines per mm gratings offer a
279 low resolution but fast measurements, and allow an overview of the Raman response over
280 the spectral range. The higher number of lines per mm offers a high resolution but longer
281 measuring times.
282 Normally the low resolution grating would be used to determine the region of interest where
283 further acquisition should be done using a higher resolution grating.

284
285 6.2.6 Sample focus
286 The quality of the spectra depends on the correct focus on the sample. In confocal mode, it
287 will allow to reduce interference from other surrounding materials. Tiny details in the sample
288 could thus be investigated using a higher magnification and confocal mode.

289
290 6.2.7 Fluorescence
291 Fluorescence is an important issue in Raman spectroscopy. It can partially or even totally
292 mask the Raman signal. Using another excitation wavelength (laser) can sometimes improve
293 the Raman response. Fluorescence problems can therefore be more easily avoided when
294 possessing multiple laser sources.
295 Another method is photo-bleaching, where the sample is irradiated by the laser for a long
296 period of time (e.g. 15 min – 1 hr) in order to reduce or suppress fluorescence. Since Raman
297 spectroscopy is already a time-consuming technique, the use of photo-bleaching may be
298 limited in routine work.
299 SERRS will also quench fluorescence and enhance the Raman signal.

300 301 6.3 Spectra redundancy and documentation

302
303 Paint samples are notoriously inhomogeneous. Analysis should be performed on a sufficient
304 number of measuring spots to account for compositional variations.

305

306 Spectra shall be stored in their raw form and unambiguously labelled, manipulated spectra
307 shall be stored with appropriate annotations. The raw file format of the instrument
308 manufacturer usually saves instrumental conditions and spectral corrections. It is
309 recommended to store the raw data in this format.

310
311

312 **7 INTERPRETATION OF RESULTS**

313
314

7.1 Spectral comparison

315
316
317
318

In a first step the homogeneity of the spectra of the questioned and the control sample shall be assessed. Big variations have to be included in the interpretation of the results.

319
320
321
322
323

For comparison purposes, baseline corrected spectra are superimposed to check for similarities and differences. If some spectra are especially noisy, it may be difficult to observe them in all detail. Smoothing is possible but must be used with care as it may affect the presence of minor bands.

324
325
326
327
328

Given spectra of good quality, the patterns of major and minor peaks, their relative positions and intensities, of the control sample are compared to the recovered paint spectra. If the spectra are comparable then the samples are said to be indistinguishable. If peaks differ substantially in position or peaks are missing in one of the compared samples, then the samples do not match.

329
330
331
332
333
334
335
336
337
338

The Raman analysis is very local, especially when using higher magnification to focus on tiny details, and paint samples are known to be inhomogeneous. In addition the Raman signal is rather proportional than quantitative and resonance effects may significantly increase the signal of the pigment compared to other components such as fillers. The relative intensities may then vary when detecting multiple components in the same spectrum. The Raman spectra from different measuring spots may consequently slightly differ and it is sometimes desirable to compare several Raman spectra together to obtain an overview of the peaks present or absent from a sample.

339
340

7.2 Spectral interpretation using Raman libraries

341
342
343
344

Raman spectra contain information from the paint binder, its pigments and fillers. The dominant bands originate from pigments and fillers. When a mixture is present, one of the constituents is often dominant and can hide contributions of the others.

345
346
347
348
349
350

The number of pigments is limited to a few hundred. The ENFSI EPG has created a Raman database using 458, 514, 633 and 785 nm lasers. This database contains almost all pigments available on the market. Each pigment has a unique Raman spectrum and its identification when used alone is relatively straightforward. When a mixture of pigments with a Raman contribution is present, identification becomes more complicated and using several laser sources may be helpful to detect more than one pigment.

351
352
353

7.3 Use of KnowItAll software

354
355

Raman spectra of pigments and fillers can be compared to the EPG Raman library of pigments by using the KnowItAll software. Basic recommendations for its optimal use are:

356
357 7.3.1. Data correction
358 Data correction steps such as normalization and baseline correction may be done temporarily
359 before the library search. It is preferable to use the same conditions for all spectra to be
360 compared. A strong baseline correction is recommended before searching in libraries.
361

362 7.3.2. Search algorithm
363 Routines for spectra, peak and property search are available and can be used separately or
364 in combination. A non-derivative search algorithm (Euclidian distance or Correlation) without
365 prompting the automated baseline correction is sufficient.
366
367

368 **8. REFERENCES**

- 369
- [1] ILAC G19:08/2014, Modules in a Forensic Science Process
 - [2] ISO/IEC 9000:2015, Quality management systems – Fundamentals and vocabulary.
 - [3] ISO/IEC 17020:2012, Conformity assessment – Requirements for the operation of various types of bodies performing inspection
 - [4] ISO/IEC 17025:2017, General requirements for the competence of testing and calibration laboratories
 - [5] ASTM E1492-11 (2017), Standard Practice for receiving, documenting, storing and retrieving evidence in a forensic science laboratory
 - [6] ASTM E1840-96 (2014) Standard Guide for Raman shift standards for spectrometer calibration

370 **General literature**

371
372
373 Dupuie J.L., Weber W.H., Scholl D.J., Gerlock J.L., Clearcoat analysis in isolated and multi-
374 layer paint systems by confocal Raman microscopy, *Polymer Degradation and Stability* (1997)
375 57: 339-348.

376
377 Massonnet G., Stoecklein W., Identification of organic pigments in coatings: applications to
378 red automotive topcoats, Part III. Raman spectroscopy (NIR FT-Raman), Proc. 1st Meeting of
379 the European Academy of Forensic Sciences, Lausanne, Switzerland (1997)

380
381 McCreery, R. L. 2000. Raman spectroscopy for chemical analysis. Chichester, UK: Wiley-
382 Interscience.

383
384 White, P. C. 2000. SERRS Spectroscopy, a new technique for Forensic Science? *Science &*
385 *Justice* 40: 113-119.

386
387 Suzuki, E. and Carrabba, M., 2001, "In Situ Identification and Analysis of Automotive Paint
388 Pigments Using Line Segment Excitation Raman Spectroscopy: I. Inorganic Topcoat
389 Pigments," *Journal of Forensic Sciences*, Vol. 46, No. 5, 0, pp. 1053-1069,

390
391 De Gelder J., Vandenabeele P., Govaert F., Moens L., Forensic analysis of automotive paints
392 by Raman spectroscopy, *J. Raman Spectrosc.* (2005) 36: 1059-1067.
393

394 Muehlethaler C., Massonnet G., Esseiva P., The application of chemometrics on Infrared and
395 Raman spectra as a tool for the forensic analysis of paints, *For. Sci. Intern.* (2011) 209: 173-
396 182.
397
398 Zieba-Palus J., Michalska A., Photobleaching as a useful technique in reducing of
399 fluorescence in Raman spectra of blue automotive paint samples, *Vibrational Spectrosc.*
400 (2014) 74: 6-12.
401
402 Trzcinska B., Kowalski R., Zieba-Palus J., Comparison of pigment content of paint samples
403 using spectrometric methods, *Spectrochim. Acta A* (2014) 130: 534-538.
404
405 Muehlethaler C., Massonnet G., Buzzini P., Influence of the shaking time on the forensic
406 analysis of FTIR and Raman spectra of spray cans, *Forens. Sci. Intern.* (2014) 237: 78-85.
407
408 Lambert D., Muehlethaler C., Esseiva P., Massonnet G., Combining spectroscopic data in the
409 forensic analysis of paint: Application of a multiblock technique as chemometric tool, *For. Sci.*
410 *Intern.* (2016) 263: 39-47.
411
412 Michalska A., Martyna A., Zadora G., Investigation of various factors influencing Raman
413 spectra interpretation with the use of likelihood ratio approach, *Forens. Sci. Intern.* (2018) 282:
414 60-73.
415
416 Muehlethaler C., Cheng Y.P., Islam S.K., Lombardi J.R., Contribution of Raman and Surface
417 Enhanced Raman Spectroscopy (SERS) to the analysis of vehicle headlights: Dye(s)
418 characterization, *Forens. Sci. Intern.* (2018) 287: 98-107.
419
420 X. He, J. Wang, C. Gao, Y. Liu, Z. Li, N. Li, J. Xia, Differentiation of white architectural paints
421 by Microscopic Laser Raman spectroscopy and Chemometrics, *Spectrochim. Acta A:*
422 *Molecular and Biomolecular Spectroscopy* (2020)
423
424 Zieba-Palus J., Examination of the variation of chemical composition and structure of paint
425 within a car body by FT-IR and Raman spectroscopies, *J. Molec. Struct.* (2020) 1219: 128558
426
427 Martyna A., Menzyk A., Damin A., Michalska A., Martra G., Alladio E., Zadora G., Improving
428 discrimination of Raman spectra by optimizing preprocessing strategies on the basis of the
429 ability to refine the relationship between variance components, *Chemometr. Intell. Lab.*
430 *Systems* (2020) 202: 104029.
431
432 Bosi A., Ciccola A., Serafina I., Guiso M., Ripanti F., Postorina p., Curini R., Bianco A., Street
433 art graffiti: discovering their composition and alteration by FTIR and micro-Raman
434 spectroscopy, *Spectrochim. Acta A: Molecular and Biomolecular Spectroscopy* (2020) 225:
435 117474.
436
437

9 AMENDMENTS AGAINST PREVIOUS VERSION

New document

###